

Upload form for Libraries/Methods

Name	HSLExtensions v1.1.0	Creation date	01.02.2013
Name of designer	Oliver Funke	E-Mail address	Ofunke@hamilton.ch
Purpose	Framework with basic libraries.		
Functions	<p>The following libraries exist: Array, Array2D, Array3D, Core, Dictionary, Directory, Drive, File, Framework, Pipetting, Sequence, SequenceArray2D, String, Venus, Windows.</p> <p>New Features:</p> <ul style="list-style-type: none"> • Validation status of all libraries is set. Appropriate test plan and test result documentation added. • Dictionary::GetKeys • Dictionary::GetItems • Sequence::Randomize • Sequence::RemovePositions • Sequence::RemovePositionsByProperty • Sequence::ReplaceLabwareIdentifier • SequenceArray2D::AddRow • SequenceArray2D::Create • SequenceArray2D::GetColumn • SequenceArray2D::GetRow • SequenceArray2D::GetSize • SequenceArray2D::GetSequence • SequenceArray2D::SetColumn • SequenceArray2D::SetRow • SequenceArray2D::SetSequence 		
Search strings <i>(only necessary for libraries)</i>			
Parameters to pass <i>(only necessary for libraries)</i>			
Decklayout/System screenshot <i>(only necessary for methods)</i>			
Instrument configuration <i>(please double click the boxes)</i>			
<input type="checkbox"/> ML StarLet <input type="checkbox"/> ML Star <input type="checkbox"/> ML StarPlus <input type="checkbox"/> Autoload <input type="checkbox"/> iSwap <input type="checkbox"/> CO-RE Grip 1000ul <input type="checkbox"/> CO-RE Grip 5ml <input type="checkbox"/> CO-RE 96 Head <input type="checkbox"/> CO-RE 384 Head <input type="checkbox"/> Nano Pipettor <input type="checkbox"/> Wash Station 1 <input type="checkbox"/> Wash Station 2 <input type="checkbox"/> TCC 1 <input type="checkbox"/> TCC 2 <input type="checkbox"/> Tube Gripper <input type="checkbox"/> Camera Channel Pump Stations: HW <input type="checkbox"/> HV <input type="checkbox"/> HU <input type="checkbox"/> # of Channels 1000ul: # of Channels 5ml: # of S-Tube Capper:			
Software & Operating System <i>(please double click the boxes)</i>			
Hamilton software	<input type="checkbox"/> V3.2.1 <input type="checkbox"/> V3.2.3 <input type="checkbox"/> V4.0.0 <input type="checkbox"/> V4.0.1 <input type="checkbox"/> V4.1.0 <input type="checkbox"/> V4.2 IVD <input type="checkbox"/> VENUS one <input type="checkbox"/> VENUS two		
Additional	<input type="checkbox"/> Scheduler <input type="checkbox"/> TADM <input type="checkbox"/> Database plus <input type="checkbox"/> SWAP		
Hamilton service pack	<input type="checkbox"/> Service Pack 1 <input type="checkbox"/> Service Pack 2 <input type="checkbox"/> Service Pack 3		
Firmware version	<i>(please type the firmware versions from tracefile)</i>		
Operating system	<input type="checkbox"/> Win2000 <input type="checkbox"/> SP1 <input type="checkbox"/> SP2 <input type="checkbox"/> SP3 <input type="checkbox"/> SP4 <input checked="" type="checkbox"/> Win XP <input type="checkbox"/> SP1 <input type="checkbox"/> SP2 <input type="checkbox"/> SP3 <input checked="" type="checkbox"/> Win VISTA <input type="checkbox"/> SP1 <input type="checkbox"/> SP2 <input checked="" type="checkbox"/> Other: Windows 7		
Third party components <i>(please list all used third party components as exact as possible)</i>			
Driver/Library list <i>(please list all drivers and libraries necessary to run the protocol)</i>			
ASW Standard Libraries			
How to install/How to uninstall			
Install: execute " Hamilton.HSLExtensions.Setup_v1.1.0.msi" Uninstall:execute " Hamilton.HSLExtensions.Setup_v1.1.0.msi" OR remove " HSLExtensions" with the windows "Add/Remove Programs" control panel			