

## Release notes – VENUS FIVE

### ML STAR:

The version 4.6.0.5313 will be released as "ML STAR VENUS five Software" (Biotech) and includes the following changes/new features:

Function	Description / Changes / New features
<b>Phoenix</b>	<ul style="list-style-type: none"><li>• Use of Phoenix 4.6.0.8060</li></ul>
<b>General</b>	<ul style="list-style-type: none"><li>• Support of Windows 10 64-bit operating system.</li></ul>
<b>ML STAR Smart Steps</b>	<ul style="list-style-type: none"><li>• Support of instrument status light</li></ul>
<b>ML STAR Steps</b>	<ul style="list-style-type: none"><li>• Support of instrument status light</li><li>• New single steps Set Status Light and Status Light Progress</li><li>• New single steps for VacuTube Channel (new Labware Handling Channel type)</li></ul>
<b>ML STAR configuration</b>	<ul style="list-style-type: none"><li>• Barcode settings: add 2D barcode types (default all disabled)</li><li>• Instrument configuration: Status Light Progress (default Off)</li><li>• Simulator configuration:<ul style="list-style-type: none"><li>- Autoload: configurable Not available/ 1D Reader / 2D Reader (default 1D Reader)</li><li>- Labware Handling Channel type: add VacuTube Channel</li></ul></li><li>• Status Light (default Available)</li></ul>
<b>Maintenance</b>	<ul style="list-style-type: none"><li>• Support of instrument status light</li></ul>
<b>Bugfixes</b>	<ul style="list-style-type: none"><li>• Bug 5032 (SCR #1667): 96er TADM Head - Tips not ejected in Clean-Up Routine</li><li>• Bug 6127 (SCR #1657): 5ml Tips not ejected in Clean-up routine</li><li>• Bug 6130 (SCR #1660): Smart Step Aliquote - Wrong liquid following in connected container</li><li>• Bug 6134 (SCR #1669): CO-RE 96 Head Move (Single Step) - Parameter out of range Error</li><li>• Bug 6136 (SCR #1671): Wrong Dispense height in Exclude Error handling after TADM clot error</li></ul>

**Phoenix:**

Version 4.6.0.8060 has following changes / new features:

Function	Description / Changes / New features
<b>General</b>	<ul style="list-style-type: none"><li>• Support of Windows 10 64-bit operating system.</li></ul>
<b>Editors</b>	<p>Hamilton Method Editor:</p> <ul style="list-style-type: none"><li>• New Quick Start dialog when method editor is launched (also accessible anytime from the File Menu).</li></ul> <p>System Deck Editor:</p> <ul style="list-style-type: none"><li>• Changed layout of system deck editor panels</li><li>• Integrated sequence positions grid and sequence sort options into sequence panel (as replacement for previously advanced dialog "Edit Positions")</li><li>• New Function to delete duplicate sequence positions</li><li>• Views not related to system deck editor are automatically hided whenever the system deck editor gets the focus</li><li>• Changed visualization of Labware search results</li></ul> <p>Labware Editor:</p> <ul style="list-style-type: none"><li>• New labware editor with support for Multiflex Carrier definition</li><li>• Configurable to use either the new Graphical Labware Editor (default) or the Legacy Labware Editor</li></ul>
<b>General Steps</b>	<ul style="list-style-type: none"><li>• File Open step: support of Microsoft Excel xlsx format</li></ul>
<b>Power Steps</b>	<p>The following Power Steps are available:</p> <ul style="list-style-type: none"><li>• Transfer Samples</li><li>• Add Reagent</li><li>• Serial Dilution</li><li>• Replicates</li><li>• Hit Picking</li></ul>
<b>Custom Dialog Step</b>	<ul style="list-style-type: none"><li>• New Template 'Hamilton'</li></ul>
<b>Bugfixes</b>	<ul style="list-style-type: none"><li>• Bug 9838: Scheduler - Max planned duration ignored</li><li>• Bug 6138 (SCR 1673): users.cfg not signed error</li><li>• Bug 6132 (SCR 1662): Memory allocation on Maintenance and Verification API</li></ul>
<b>Pre-Requirement Installer</b>	<ul style="list-style-type: none"><li>• Activation of .NET 3.5 feature (internet access required) and check for .NET 4.6.2 availability</li><li>• Microsoft® SQL Server 2019 Express Edition CU13 (optional)</li></ul>