

Project-Name :	Phoenix	HAMILTON Company	2017-03-31	Page 32 of 34
Doc.-Title :	Releases of BDZ Phoenix Products	Doc.-No. :	E2830073m	Doc.-Version : 4.5

### 3.17 Version 4.5.0.5217

The version 4.5.0.5217 was released as two different products: "ML STAR VENUS four Software" (Biotech) and "ML STAR IVD Software". The version includes the following changes/new features:

Function	Description / Changes / New features
Phoenix	<ul style="list-style-type: none"> <li>• Use of Phoenix 4.5.0.7977</li> </ul>
Firmware	<ul style="list-style-type: none"> <li>• See [3] and [4]</li> </ul>
General	<ul style="list-style-type: none"> <li>• Support of Windows 7 and Windows 10 operating systems.</li> </ul>
ML STAR Smart Steps	<ul style="list-style-type: none"> <li>• Option for customizable error handling</li> <li>• Option for sequence merging within Pipette Smart Steps</li> </ul>
ML STAR Easy Steps	<ul style="list-style-type: none"> <li>• Implemented SCR 1593: Improper Aspiration (Error Recovery - Waste)</li> </ul>
ML STAR Steps	<ul style="list-style-type: none"> <li>• AutoLys Steps error recovery with option to release gripped AutoLys Tubes (user intervention required).</li> <li>• Implemented SCR 1598: Tube Twister Channel - Labware Exist Dialog (clean-up routine at method end)</li> <li>• Implemented SCR 1601: Dispense on the Fly – Unhandled exception (loop over 12 nunc 96 plates)</li> <li>• Implemented SCR 1609: Sequence index out of range error during aliquot dispense (Smart Step – Pipette Aliquot)</li> <li>• Implemented SCR 1623: AutoLys Channel unhandled error - No element found (Firmware error code 21/88)</li> <li>• Implemented SCR 1633: Tip Eject into Liquid Waste</li> <li>• Implemented SCR 1638: Eject of CO-RE Grip Tools fails with step lost</li> <li>• Implemented SCR 1641: Repeat Error Handling - Executed twice</li> </ul>
AT File Filter	<ul style="list-style-type: none"> <li>• Implemented SCR 1614: AT File Filter - Wrong update of Register.inf file</li> </ul>
HSL Libraries	<ul style="list-style-type: none"> <li>• Implemented SCR 1627: GetBarcode of HSLMStarStepReturnLib trims leading and trailing spaces</li> </ul>
Maintenance	<ul style="list-style-type: none"> <li>• Implemented SCR 1610: Punch Gripper crash after daily maintenance</li> </ul>