

Library upload form

Name of Library	FileDirectoryLib
Name of designer	Yvonne Dangendorf
Purpose of the library	This library is for a better handling with files and directories
Creation date	2009/03/09
Tested on	
Vector software version	Venus One
Service Pack of Vector software	-
Operating system	Microsoft Windows XP Professional Version 2
Service Pack(s) of OS system	Service Pack 2
Instrument	-
Configuration	-
Firmware version(s)	-
Parameters to pass	
Search strings	File, Directory
Abstract	The HSL FileDirectoryLib provides the handling with data files and directories. Functions like deleting, renaming, copying moving and getting information of files and directories and more are included.
Functions	<ul style="list-style-type: none"> - CheckIfDirectoryExists - CheckIfFileExists - CopyFile - CreateDirectory - DeleteAllFilesInDirectory - DeleteDirectory - DeleteFile - GetFilesInDirectory - GetNumberOfFilesInDirectory - GetSubdirectoriesOfDirectory - MoveDirectory - MoveFile - RenameFile
How to install	<ul style="list-style-type: none"> - Get file HamiltonOneChoiceSetup.msi - Install as usual
How to uninstall	<ul style="list-style-type: none"> - Systemconfiguration/Software - Uninstall as usual