

# Library upload form

Name of Library	HSLXLib
Name of designer	Beni Andenmatten
Purpose of the library	Handling of Excel Files
Creation date	
Tested on	
Vector software version	4.1.0.5617
Service Pack of Vector software	-
Operating system	Microsoft Windows XP Professional Version 2
Service Pack(s) of OS system	Service Pack 2
Instrument	-
Configuration	-
Firmware version(s)	-
Parameters to pass	-
Search strings	Excel, XLib
Abstract	The HSLXLib Library provides the handling with data files provided by Excel
Functions	<ul style="list-style-type: none"> <li>- ClearCell</li> <li>- ClearCells</li> <li>- Close</li> <li>- GetActiveCell</li> <li>- New</li> <li>- Open</li> <li>- ReadCell</li> <li>- ReadCell2</li> <li>- RunMacro</li> <li>- RunMacro1</li> <li>- RunMacro2</li> <li>- RunMacro3</li> <li>- SetVisible</li> <li>- WriteCell</li> <li>- WriteCell2</li> </ul>
How to install	- Copy all files into your Library folder
How to uninstall	