

# Library upload form

Name of Library	HSLMicrosoftExcel
Name of designer	Joachim Schmid
Purpose of the library	File handling with Excel files
Creation date	2004-11-15
Tested on	
Vector software version	4.2.0.6425
Service Pack of Vector software	-
Operating system	Microsoft Windows XP Professional Version 2
Service Pack(s) of OS system	Service Pack 2
Instrument	-
Configuration	-
Firmware version(s)	-
Parameters to pass	-
Search strings	Excel
Abstract	The HSLMicrosoftExcel Library provides the handling of Excel Workbooks and Sheets. It is possible to count the sheets of a workbook, to get the name of the sheets and to read out and write into specific cells (by row and column).
Functions	<ul style="list-style-type: none"> <li>- Init</li> <li>- End</li> <li>- SetCellValue</li> <li>- GetCellValue</li> <li>- GetSheetCount</li> <li>- GetSheetName</li> <li>- SelectSheet</li> <li>- SafeFile</li> <li>- OpenFile</li> </ul>
How to install	- Execute HSLExcLib.msi
How to uninstall	- Uninstall HSLExcLib.msi via the Control Panel -> Add or Remove Programs